

NEW NURSECALL SYSTEM INFO

BEDROOM EMERGENCY CALLSTATION



SF521UL Room Station



SF301B Call Cord

DO NOT UNPLUG CORD!!!
ONLY PRESS IN A MEDICAL EMERGENCY!!!

IN THE EVENT OF A **MEDICAL EMERGENCY** - PUSH THE BUTTON ON THE CORD TO CALL FOR EMS.

BATHROOM EMERGENCY CALLSTATION



SF523UL Room Station

IN THE EVENT OF A **MEDICAL EMERGENCY** - PUSH THE RED BUTTON ON THE LOWER HALF OF THE ROOM STATION TO CALL FOR EMS.

AFTER THE EMERGENCY BUTTON HAS BEEN PUSHED, THE RESET BUTTON **MUST BE PUSHED** TO CLEAR THE EMERGENCY FROM THE SYSTEM. FAILURE TO DO SO MAY RESULT IN NO EMERGENCY SIGNAL TRANSMISSION TO THE EMS DISPATCHER.



WHAT HAPPENS AFTER AN EMERGENCY STATION BUTTON IS PUSHED?

- THE NURSE CALL SYSTEM SENDS A SIGNAL VIA PHONE LINES TO A CENTRAL MONITORING STATION. THE CENTRAL MONITORING STATION THEN CONTACTS THE LOCAL FIRE/EMS FOR A MEDICAL DISPATCH.

WHAT INFORMATION IS TRANSMITTED WHEN AN EMERGENCY BUTTON IS PUSHED?

- THE ADDRESS OF THE BUILDING IS TRANSMITTED.
- THE UNIT NUMBER IS TRANSMITTED.
- THE NAME AND PHONE NUMBER OF THE PRIMARY PERSON ON RECORD IS TRANSMITTED.
- EMS WILL RESPOND WITH THE ABILITY TO ACCESS THE BUILDING THROUGH THE MAIN ENTRY VIA A LOCK-BOX THAT YOUR PROPERTY MANAGER HAS COORDINATED WITH EMS SERVICES.